

MICRO Q TECHNOLOGIES

BOVINE

PORTABLE TRANSPORTABLE INCUBATOR

PORTABLE TRANSPORTABLE INCUBATOR

PLATFORMS FEATURE

- Active Ultra-Precision Heating
- Temperature control range 20° C to 45° C
- Ambient Temperature range -10° C to 45° C
- Microprocessor Controlled Heat Pump
- Temperature Control Resolution = 0.0625° C



- Temperature LCD Display Resolution = 0.1° C
- Rechargeable Battery with Integrated Power Supply
- Temperature Holding Duration 72+ hours
- Vibration / Shock Resistant Exterior
- Extremely Durable Ballistic Nylon Skin

PORTABLE TRANSPORTABLE INCUBATOR

Micro Q iQ1 T provides ultra-precise temperature control of oocytes and embryos during transport. The incubator is extremely compact, light and very flexible with customization.



TRANSPORTABLE EMBRYO INCUBATOR DIMENSIONS

- Length: 10.0 inches
- Width: 9.50 inches
- Height: 10.0 inches
- Weight: 7 pounds



LOCATION

14362 North Frank Lloyd Wright Blvd Suite 1210
Scottsdale, Arizona 85260 USA



CALL US

Phone: 480.380.9669



EMAIL US

sales@microq.com